


## Spare Parts List Incubator 8000 IC, US version



This spare parts list is valid in the USA and Canada for part no.:

Part No.	Description
<b>2M20550</b>	INCUBATOR 8000 IC USA
<b>2M21250</b>	INCUBATOR8000 IC USA SCAN.HOOD
<b>2M21966</b>	BASIC UNIT INC.8000IC V - 120V
<b>2M21973</b>	INCU.8000 IC USA VLH 120V SDW
<b>FR00054</b>	VARI.-SYSTEM INCUBATOR 8000 IC


## Summary of pictures

Picture	Description
2	BASIC UNIT INCUBATOR 8000 IC V-- 120V
3	BASIC UNIT INCUBATOR 8000 IC V-- 120V
4	CARRIAGE WITH LIFTING COLUMN 115V
5	PEDAL, CPL.
6	PLATE, CPL.
7	FAN WHEEL, CPL.
8-11	HOUSING NC, IC, SC, CPL.
12	HINGE BRACKET FOR HOOD
13	TAPPET,CPL., FOR STANDARD HOOD
14	CATCH, FOR SCAND. HOOD
15	O2 SUPPLY INCUBATOR 8000
16	MODIF.KIT PRESSURE REGULATOR
18	SWIVEL CUPBOARD, CPL.
19	MOUNTG.DEVICES F.CUPB. INCUBATOR
22-23	ELECTRIC MAINS CONNECTION IC
24	AGGREGATE
25	AGGERATE SHEET METAL WITH PCB, CPL.
26	MOTOR, CPL.
27	HUMIDIFIER, CPL.
28	HANDLE, CPL.
29	INCLINATION, RIGHT/LEFT
30-31	MODIFICATION SET H2O SUPPLY / INCUBATOR 8000 SC TO IC
32	WATER CONTAINER SET (FROM INCUBATOR 8000 SC)
33	HOOD INC. 8000 IC, CPL., STANDARD
34	HINGE,CPL., SENSOR
35	FLAP, CPL., FOR STANDARD HOOD
36	SET OF SWIVELLING WINDOWS

Picture	Description
37	HOOD INCU. 8000 IC SCAND.
38	FLAP, BIG (CPL.), FOR SCAND. HOOD
39-40	FLAP, SMALL (CPL.), FOR SCAND. HOOD
44	SENSOR, COMPL., INCUBATOR 8000 IC
45	BORDERING, CPL.
46-52	PLUG-IN UNIT 8000 IC
53	PCB DISPLAY AIR
54	PCB POWER PACK
55	PCB CPU INCUBATOR
56	PCB ANALOG
57	PCB DISPLAY / DISPLAY SKIN
58	INSTALLATION SET RS 232
60	THERMOMONITORING, CONNECTION FOR PROBES
62	SKINTEMPERATURE SENSORS
63	MODIFICATION SET, INCUBATOR 8000 IC, US LANGUAGE
65	VENTILATION HOSE HOLDER
66	MODIFICATION SET INCUBATOR IC, USA
69	MODIFICATION KIT CORNER PIECE ACCORDING TO TSB # 4

## BASIC UNIT INCUBATOR 8000 IC V-- 120V


Picture 2



Datei/File: 6141232.002

**BASIC UNIT INCUBATOR 8000 IC V-- 120V**

**Picture 3**




Datei/File: 6141232.003

Item No.	Description	Order-Code
1	MAINS CABLE UL	2M20453
2	CARRIAGE WITH LIFTING COLUMN 120V, <a href="#">SEE PICTURE 4</a>	
3	FUSE LINK DIN 41662 T6,3A	1815172
4	AGGREGATE, CPL., <a href="#">SEE PICTURE 24</a>	2M20615
6	SUPPORT, CPL.	2M19654
8-10	MODIFICATION KIT CORNER PIECE, <a href="#">SEE PICTURE 69</a>	2M20722
11	HANDLE,CPL., <a href="#">SEE PICTURE 28</a>	2M20391
13	HANDWHEEL, WHITE	2M21202
13	HANDWHEEL, ORANGE	2M19647
14	SUPPORT, LEFT	2M20692
15	SCREW DWN 562-4X20-C-ST/A2P	1338382
16	WASHER B 4,3 DIN 9021-A4	1331914
17	SUPPORT, RIGHT	2M20660
18	INCLINATION, RIGHT, <a href="#">SEE PICTURE 29</a>	2M19712
20	SPLIT WASHER B6 DIN127-STAINL.	1333941
22	INCLINATION, LEFT, <a href="#">SEE PICTURE 29</a>	2M19650
23	FAN WHEEL, CPL., <a href="#">SEE PICTURE 7</a>	2M20205
26-27	LYING SURFACE	2M20888
28	MATTRESS	2M20907
29	INTERMEDIATE ELEMENT	2M19537
30	TUB	2M19334
31	GASKET	2M19595
32	HOUSING NC,SC,IC, CPL., <a href="#">SEE PICTURE 8</a>	
33	O2 SUPPLY, <a href="#">SEE PICTURE 15</a> AND <a href="#">SEE PICTURE 16</a>	2M20666
34	BED EXTENSION	2M21225
35	BORDERING,CPL., <a href="#">SEE PICTURE 45</a>	
38	ELECTRIC MAINS IC, <a href="#">SEE PICTURE 22</a>	
39	CABLE CLIP (2,4X92)	8712007
40	BLUE NOSE CAT, without illustration	2M21420

Item No.	Description	Order-Code
42	MATRESS-CLOTH, without illustration	2M21272
45	SCREW DWN 562-4X20-C-ST/A2P	1338382
50-55	<a href="#">SEE PICTURE 69</a> and TSB # 4	

## CARRIAGE WITH LIFTING COLUMN 230V


Picture 4



Datei/File: 6141232.004



Item No.	Description	Order-Code
1-5	HEIGHT ADJUSTMENT SC/IC 100/127V 50/60HZ	2M20244
1	ALTIT. ADJUSTMENT SC/IC 100/127V 50/60HZ	2M20296
2	PEDAL, CPL., <a href="#">SEE PICTURE 5</a>	2M20290
4a	FIXING GUIDE ROLL 125 MM	2M20794
5a	GUIDE ROLL 125 MM	2M20260
4b	FIXING GUIDE ROLL 100 MM	2M20792
5b	GUIDE ROLL 100 MM	2M20250
6-21	CARRIAGE, CPL.	2M20240
11	WASHER B 6,4 DIN 125-A4	1330705
12	SPLIT WASHER B6 DIN127-STAINL.	1333941
13	CH. HEAD SCREW M4X8 DIN912	1263056
14	SPLIT WASHER B4 DIN 127	1331345
15	FLAT PLUG DIN 46244-A6,3	1814079
18	CAP, WHITE	2M21223
18	CAP, ORANGE	2M20348
19	TUBE, BACK	2M20232
20	TUBE, FRONT	2M20234


**PEDAL, CPL.****Picture 5**

Datei/File: 6141232.005

Item No.	Description	Order-Code
1-7	PEDAL, CPL.	2M20290
1	PEDAL	2M20218
2	PLATE, CPL.	2M20370
4	WASHER	G80073
5	SELF-LOCKING NUT M6 DIN 980-A2	1333062


## PLATE, CPL.

Picture 6




Datei/File: 6141232.006

Item No.	Description	Order-Code
1-7	PLATE, CPL.	2M20370
2	MICROSWITCH	1831666
5	FLAT PLUG DIN 46244-A6,3	1814079
6	LOCK WASHER	1293680
7	SCREW AM 4X10 DIN84-A4/051	1329049

**FAN WHEEL, CPL.****Picture 7**

Datei/File: 6141232.007


Item No.	Description	Order-Code
1-3	FAN WHEEL, CPL.	2M20205
1	CONE	2M20542
2	FAN WHEEL	2M20203
4	TOROIDAL SEALING RING	M19241

**HOUSING NC, IC, SC, CPL.****Picture 8**

Datei/File: 6141232.008

HOUSING NC, IC, SC, CPL.


Picture 9



Datei/File: 6141232.009

**HOUSING NC, IC, SC, CPL.**


### Bild/Picture 10



Datei/File: 6141232.010

HOUSING NC, IC, SC, CPL.


Picture 11



Datei/File: 6141232.011




Item No.	Description	Order-Code
2	SCREW DWN 562-4X20-C-ST/A2P	1338382
4	SCREW BZ 2,9X19 DIN 7971-A2	1337610
9	CONNECTOR	2M20332
15	LOCK	2M20257
16	LABEL	2M20629
17	HOUSING	2M20299
18	GASKET	2M20199
20	HINGE	2M19730
23	CABLE CLIP (2,4X92)	8712007
24	VALVE SEAT	RM08303
25	CONTROL VALVE PLATE	RM08304
26	TOROIDAL SEALING RING	D11137
28	FILTER	8402926
29	HINGE BRACKET, CPL., <a href="#">SEE PICTURE 12</a>	
30	TAPPET,CPL., <a href="#">SEE PICTURE 13</a>	2M19543
33	FAN	2M20021
34	EARTH CONDUCTOR TERMINAL	8300390

**HINGE BRACKET, CPL.****Picture 12**


Datei/File: 6141232.012

Item No.	Description	Order-Code
1	HINGE BRACKET	2M19542
2	SCREW	2M19529
3	LOCK WASHER 7 DIN 6799	1309773

**TAPPET, CPL. FOR STANDARD HOOD****Picture 13**


Datei/File: 6141232.013

Item No.	Description	Order-Code
1-7	TAPPET,CPL.	2M19543
3	LOCK WASHER 12 DIN 6799	1291467
4	THRUST WASHER 4 DIN 6799	1328999
5	GUIDE	2M19531
6	SPRING	T50137

**CATCH FOR SCAND. HOOD****Picture 14**

Datei/File: 6141232.033

Item No.	Description	Order-Code
1-5	SPRING-LOADED CATCH	2M20852
1	HOSE 5X2SI 60SHA NF M25779	1203606
2	BOLT	2M20853
3	LOCK WASHER 12 DIN 6799	1291467
4	GUIDE	2M19531
5	THRUST WASHER 4 DIN 6799	1328999


**O2 SUPPLY****Picture 15**

Datei/File: 6141232.014

Item No.	Description	Order-Code
1-16	O2 SUPPLY	2M20666
4	2/2-PORT DISTRIBUTING VALVE	2M19958
7	WASHER B 4,3 DIN 9021-A4	1331914
11	WASHER A 4,3 DIN 125-A2	1327542
12	SCREW DWN 562-4X20-C-ST/A2P	1338382
13	MODIF. KIT PRESSURE REGULATOR, <a href="#">SEE PICTURE 16</a>	
14	ADHESIVE LABEL	8303296
16	CABLE CLIP (2,4X92) without illustration	8712007


## MODIF.KIT PRESSURE REGULATOR

Picture 16



Datei/File: 6141232.015


Item No.	Description	Order-Code
1	PRESSURE REGULATOR	8402745
3	TUBE 6X1 DIN73378-PA 12 W NF.	1202618
4	THRUST COLLAR 6, WHITE	M31603
5	ANGLE CONNECTION	M30961
8	FILTER	2M19622
9	SCREW DWN 562-4X20-C-ST/A2P	1338382
	O2 CONNECTION (USA)	8409401

**SWIVEL CUPBOARD, CPL.**
**Picture 18**


Datei/File: 6141232.016

Item No.	Description	Order-Code
1-6	SWIVEL CUPBOARD, CPL.	2M20638
1	SWIVELLING DRAWER	2M20627
2	INSERT	2M20642
3	SWIVELLING PLATE	2M20628
4	TUBE	2M20639
5	CATCH	2M20632
6	SPRING	2M20645




**MOUNTG. DEVICES F. CUPB. INCUBATOR****Picture 19**

Datei/File: 6141232.017

Item No.	Description	Order-Code
1-6	MOUNTG. DEVICES F. CUPB. INCUBATOR	2M20656
1	NUT	1267582
3	SCREW DIN 912-M8X390+-2 8.8/155	1342096
4	DISK SPRING A25 DIN 2093	1270338
5	B8,4 DIN 9021-A4	1335367


## Picture 22



26/87

## ELECTRIC MAINS IC

Picture 23




Datei/File: 6141232.059

Item No.	Description	Order-Code
6	CABLE CLIP (2,4X92)	8712007
10	SCREW AM3X10 DIN 85-A2	1333984
11	HEXAGON NUT M3 DIN934-M A4	1334913
13	SUPPORT WITH CAP	1832662
17	SPLIT WASHER B4 DIN 127	1331345
18	WASHER A 4,3 DIN 125-A2	1327542
19	SCREW DWN 562-4X20-C-ST/A2P	1338382
41	FUSE LINK DIN 41662 T6,3A	1815172
26	SHEET METAL	2M20461
29	FLAT HEAD SCREW AM 2X4 DIN84	1287664
34	HEXAGON NUT M4 DIN 934-A4/051	1328956
36	SCREW DWN 561-4X20-C-A2-H	1338072
37	FLAT PLUG DIN 46244-A6,3	1814079
38	BOLT	2M20699
39	SPLIT WASHER B3 DIN127-X12CRNI	1309749
41	FUSE LINK DIN 41662 T6,3A	1815172
42	MAINS CABLE BUSHING	2M18642

AGGREGATE, CPL.

Picture 24




Datei/File: 6141232.019

Item No.	Description	Order-Code
1-49	AGGREGATE, CPL. with humidifier, heating cartridges, fuses, fan wheel	2M20615
6	SCREW AM4X6 DIN 84-A4	1330829
7	SPLIT WASHER B4 DIN 127	1331345
8	CABLE CLIP (2,4X92)	8712007
10	WASHER B 4,3 DIN 9021-A4	1331914
11	HEXAGON NUT M4 DIN 934-A4/051	1328956
12	TEMPERATURE FUSE 315 DEGREE	1823809
13	SILICONE HOSE 6X2,5 NF M29907	1197851
15	DISC B 5,3 DIN 9021-A4	1329472
16	SCREW AM 4X12 DIN 84-A2	1287710
20	GASKET	2M19637
22	HEXAGON NUT M6 DIN 934-M A4	1267515
24	MOTOR,CPL., <a href="#">SEE PICTURE 26</a>	
25	SPLIT WASHER B3 DIN127-X12CRNI	1309749
26	HEXAGON NUT M3 DIN934-M A4	1334913
27	WASHER	1330233
29	PACKING RING	2M19638
31	COIL	2M18839
33	HUMIDIFIER, <a href="#">SEE PICTURE 27</a>	
34	ADHESIVE LABEL, without illustration	8303296
37	FAN WHEEL, CPL., <a href="#">SEE PICTURE 7</a>	
38	HEATING CARTRIDGE (120/127V) FOR HUMIDIFIER, <a href="#">SEE PICTURE 27</a>	
39	HEATING CARTRIDGE 120/127 V	2M20331
40 - 42	Fuse T 6,3 A DIN 41662	1815172
40 - 42	Fuse T 1A DIN41662	1815083
40 - 42	Fuse T 80 mA DIN 41662	1814974
40 - 42	Fuse T 200 mA DIN 41662	1815016
40 - 42	Fuse T 500 mA DIN 41662	1815059
40 - 42	Fuse T 2,5 A DIN 41662	1815121

Item No.	Description	Order-Code
40 - 42	Fuse T 1,6 A DIN 41662	1815105
40 - 42	Fuse T 160 mA DIN 41662	1815008
46	CONNECTION, without illustration	2M20464
50	CABLE HARNESS (PLUG-IN UNIT - AGGREGATE) without illustration	8200759

## AGGREGATE SHEET METAL, CPL.

Picture 25




Datei/File: 6141232.020



Item No.	Description	Order-Code
1	RELAY SEMICONDUCTOR	1831399
3	SPLIT WASHER B4 DIN 127	1331345
4	WASHER A 4,3 DIN 125-A2	1327542
5	EARTH CONDUCTOR TERMINAL	8300390
6	SCREW M3X6 DIN 84-A4	1330810
7	SPLIT WASHER B3 DIN127-X12CRNI	1309749
9	TRANSFORMER	8290426
10	SCREW DIN 84-AM4X50-A4	1287753
11	CABLE CLIP (2,4X92)	8712007
12	HEADER PROFILE 1-2CH15956 SW	1191640
14	SCREW	1329987
15	WASHER	1330233
17	PCB AGGREGATE	8290511
17a	RELAY 12V/8A 1XUM, (K1/K5 ON AGGREGATE PCB)	1830732
17b	RELAY 24V/6A/250V, (K2/K5 ON AGGREGATE PCB)	1821474
17c	RELAY, (K4 ON AGGREGATE PCB)	8301198
19	PCB HL-RELAY	8306151
21	B4 DIN 137-STAINLESS STEEL	1335308


**MOTOR,CPL.**

**Picture 26**




Datei/File: 6141232.021

Item No.	Description	Order-Code
1-22	MOTOR,CPL. (ORDER IN ADDITION CONE 2M20542 FOR MOTOR REPLACEMENT)	2M20350
5	MOTOR (ORDER IN ADDITION CONE 2M20542 FOR MOTOR REPLACEMENT)	2M20379
8	GASKET	2M19591
9	RUBBER SUSPENSION	2M18147
11	SPRING	E20215
12	CABLE CLIP (2,4X92)	8712007
13	NUT BM3 DIN 439-A4/051	1334743
15	SCREW DWN 561-4X20-C-A2-H	1338072
18	HEXAGON NUT M3 DIN934-M A4	1334913
23	BOTTLE OF SPECIAL OIL, without illustration	2M07839
29	SCREW AM4X6 DIN 84-A4	1330829
30	WASHER A 4,3 DIN 125-A2	1327542
31	SPLIT WASHER B4 DIN 127	1331345

**HUMIDIFIER, CPL. 120V****Picture 27**


Datei/File: 6141232.022

Item No.	Description	Order-Code
1-14, 16,18	HUMIDIFIER, CPL., 120/127V	2M20619
1	CONNECTION	2M20292
2	O-RING SEAL	2M08777
5	SCREW AM 4X12 DIN 84-A2	1287710
6	SPLIT WASHER B4 DIN 127	1331345
7	WASHER A 4,3 DIN 125-A2	1327542
10	COOKER	2M20291
11	TEMPERATURE LIMITER	2M20382
12	TEMPERATURE SENSOR, CPL.	2M20381
16	HEATING CARTRIDGE (120/127V)	2M20328
18	SCREW AM4X6 DIN 84-A4	1330829

**HANDLE, CPL.****Picture 28**

Datei/File: 6141232.023


Item No.	Description	Order-Code
1-5	HANDLE, CPL.	2M20391
2	END PIECE	2M20696
3	TUBE	2M13388
4	SCREW M 8X55 DIN 912-A2	1329693
3,5-9	HANDLE,CPL.	2M20391
7	B8,4 DIN 9021-A4	1335367
9	CLOSING PLUG	G60036

**INCLINATION, RIGHT / LEFT****Picture 29**


Datei/File: 6141232.024

Item No.	Description	Order-Code
1-19	INCLINATION, RIGHT	2M19712
1-19	INCLINATION, LEFT	2M19650
6	SPRING	2M19648
8	WASHER	CH09520
12	CHAIN	2M19660
18	RECEIVER, CPL.	2M19657



**MODI. SET H2O / INC. 8000 SC TO IC****Picture 30**

Datei/File: 6141232.025


**MODI. SET H2O / INC. 8000 SC TO IC****Picture 31**

Datei/File: 6141232.026

Item No.	Description	Order-Code
1-27	MODI. SET H2O / INC. 8000 SC TO IC	2M20820
1	SCREW DWN 561-4X20-C-A2-H, without illustration	1338072
2	BUFFER	6800478
3	BUTTON	2M19693
5	SCREW M4X12 DIN 427-A2	1332007
8	DISC B 5,3 DIN 9021-A4	1329472
9	WASHER B 4,3 DIN 9021-A4	1331914
12	GASKET (USA/CAN)	2M20635
13	KNURLED NUT	2M06004
14	BUTTON	2M12044
15	TRACTION ROPE	2M12045
16	SPRING-LOADED CATCH	2M19271
19	WATER TANK	2M19561
20	CONNECTION	2M20464
22-27	SHEET METAL	2M20461
22	CAP	2M20459
23	FLAT PLUG DIN 46244-A6,3	1814079
24	SPLIT WASHER B4 DIN 127	1331345
25	HEXAGON NUT M4 DIN 934-A4/051	1328956
26	ADH. TAPE 12X1,6 PVC-FOAM SW X	1202421
27	SHEET METAL	2M20549

WATER CONTAINER SET, INC. SC

Picture 32




Datei/File: 6141232.027

Item No.	Description	Order-Code
1-11	WATER CONTAINER SET	2M21747
1	SCREW DWN 561-4X20-C-A2-H	1338072
2-5	WATER CONTAINER,SET OF SPARES	2M20644
2	CONNECTION	2M20670
3	WATER TANK	2M20210
4	COVER	2M20211
6	CLAMP	2M20547
7	SPLIT WASHER B4 DIN 127	1331345
8	HEXAGON NUT M4 DIN 934-A4/051	1328956
9	FLAT PLUG DIN 46244-A6,3	1814079
11	ADH. TAPE 12X1,6 PVC-FOAM SW	1202421


HOOD INC. 8000 IC, CPL.

Picture 33



Datei/File: 6141232.028

Item No.	Description	Order-Code
1-32	HOOD INC. 8000 IC, CPL.	2M20510
1	HINGE, CPL., <a href="#">SEE PICTURE 34</a>	
4	SET COMBINATION PLEXIGLASS HOOD INCUBATOR IC/AC/NC (NOT FOR SCAN. HOOD)	2M22153
5	FLAP, CPL.	2M19525
6	CATCH, CPL.	2M20358
8-10	SWIVELLING WINDOWS, <a href="#">SEE PICTURE 36</a>	
11	CABLE CLIP	2M20376
13	CAP NUT	2M19484
14	HINGE	2M19730
16	POINTER BUTTON	2M20720
17	CATCH	2M19724
18	INTERIOR WALL	2M19512
19	HOSE THROUGH-BEARING	2M19511
20	STIFFENING	2M20214
21	RECEIVER FOR REMOVEABLE INNER WALLS	2M21449
22	BOLT FOR REMOVEABLE INNER WALLS	2M21448
25	LATERAL WALL	2M20508
26	STIFFENING	2M20568
31	CATCH, LEFT	2M19505
32	CATCH, RIGHT	2M19506
33	ELBOW, <a href="#">SEE PICTURE 69</a>	
	BED EXTENSION, <a href="#">SEE PICTURE 2</a>	

**HINGE, CPL.****Picture 34**


Datei/File: 6141232.029

Item No.	Description	Order-Code
1-5	HINGE,CPL.	8200596
3	8 MM III DIN 5401-X 40 CR 13	1298844



FLAP, CPL..

Picture 35




Datei/File: 6141232.030

Item No.	Description	Order-Code
1-13	FLAP, CPL.	2M19525
2-5	CATCH, CPL.	2M19899
2	CATCH	2M19482
3	HOUSING	2M19483
4	SPRING	8400496
5	PLATE	2M19481
6	CAP NUT	2M19484
	PLATE FOR CATCH 2M19899	2M19472
7	CATCH	2M19516
8	HINGE	2M19729
9	INTERIOR WALL	2M19512
10	FLAP	2M20733
11	SCREW DWN 562-4X14-F-A2-H	1342568
13	HINGE	2M20730

## SET OF SWIVELLING WINDOWS

Picture 36




Datei/File: 6141232.031

Item No.	Description	Order-Code
3,5-18	SET OF SWIVELLING WINDOWS, COMPLETE, LEFT AND RIGHT	2M19550
1-3	FLAP,LEFT (WITH HOSE BUSHING)	2M20552
2-4	FLAP,RIGHT(WITH HOSE BUSHING)	2M20545
2	TUBE BUSHING	2M20434
3	SCREW M6X6 DIN551-PA 6.6 NF	1342266
5	CATCH	2M19519
6	SCREW BZ 2,9X19 DIN 7971-A2	1337610
7	SCREW BZ 2,9X16 DIN7971-A2	1336630
8	BASE PLATE	2M19521
9	SUPPORT, TOP LEFT	2M19473
10	FLAP, LEFT	2M19477
11	FLAP, RIGHT If item 10 or 11 are ordered singly, item 3, screw 1342266, must also be ordered	2M19478
12	SUPPORT, TOP RIGHT	2M19475
13	BUFFER	2M19524
14	SUPPORT, BOTTOM RIGHT	2M19476
15	SUPPORT, BOTTOM LEFT	2M19474
16	GASKET, LEFT	2M19469
17	GASKET, RIGHT	2M19470
18	SPRING	2M20537


## HOOD INCU.8000 IC SCAND.

Picture 37



Datei/File: 6141232.032


Item No.	Description	Order-Code
1-32	HOOD INCU.8000 IC SCAND.	2M20610
1	HINGE,CPL., <a href="#">SEE PICTURE 34</a>	8200596
4	HOOD IC, SCAND.	2M20611
5	FLAP, BIG (CPL.), <a href="#">SEE PICTURE 38</a>	
7	SWIVELLING WINDOW, RIGHT, <a href="#">SEE PICTURE 36</a>	
8	CATCH, CPL.	2M20358
9	SCREW DWN 562-3X16-F-A2-H	1342630
10	SWIVELLING WINDOW, LEFT, <a href="#">SEE PICTURE 36</a>	
11	HINGE	2M19730
13	POINTER BUTTON	2M20720
14	FLAP, SMALL (CPL.), <a href="#">SEE PICTURE 39</a>	
15	CATCH	2M19724
16	INTERIOR WALL	2M19512
17	SET OF SWIVELLING WINDOWS, <a href="#">SEE PICTURE 36</a>	2M19550
18	CABLE CLIP	2M20376
19	CAP NUT	2M19484
21	HOSE THROUGH-BEARING	2M19511
22	STIFFENING	2M20214
25	RECEIVER FOR REMOVEABLE INNER WALLs	2M21449
26	BOLT FOR REMOVEABLE INNER WALLS	2M21448
27	LATERAL WALL	2M20508
30	SPRING-LOADED CATCH	2M20852
	BED EXTENSION, <a href="#">SEE PICTURE 2</a>	

**FLAP, BIG (CPL.)**
**Picture 38**


Datei/File: 6141232.034

Item No.	Description	Order-Code
1-14	FLAP, BIG (CPL.)	2M20607
1	BUFFER	2M19524
2	FLAP, BIG	2M20604
3	CATCH	2M19516
5	CAP NUT	2M19484
6-9	CATCH,CPL.	2M19899
6	PLATE	2M19481
7	SPRING	8400496
8	HOUSING	2M19483
9	CATCH	2M19482
10	HINGE	2M19729
11	SET OF SWIVELLING WINDOWS, <a href="#">SEE PICTURE 36</a>	2M19550




**FLAP, SMALL (CPL.)****Picture 39**

Datei/File: 6141232.035

FLAP, SMALL (CPL.)

Picture 40




Datei/File: 6141232.060

Item No.	Description	Order-Code
1-9	FLAP, SMALL (CPL.)	2M20612
7,10-15	FLAP, SMALL (CPL.)	2M20612
2	FLAP, SMALL	2M21471
3	CATCH	2M19516
7	HINGE	2M19729
12	BZ 2,9X13 DIN 7971-A2	1334948
14	O-RING SEAL	R50543
15	AXLE	2M19848


SENSOR, COMPL. INCU.8000 IC

Picture 44




Datei/File: 6141232.039


Item No.	Description	Order-Code
1-20	SENSOR, COMPL. INCU.8000 IC, for Thermomonitoring, for SW 21.n only, also known as 2M21680 by mistake	2M21688
1-20	SERVICE EXCHANGE SENSOR (2M21688)	SE2M21688
1-21	SENSOR CONNECTION (L/H,F,O2), old, for SW 20.n only	8290380
1-21	SERVICE EXCHANGE SENSOR (8290380)	SE8290380
2	HOSE 4X1,5-SI 50 SH A NF	1190520
3	SENSOR HOUSING	8200605
4	CABLE HARNESS, SENSOR PLUG	8290382
6	CAP, CPL.	8200618
7b	SOCKET,CPL. (for 8290380 / SW 20.n only)	8200608
9	PRINTED CIRCUIT BOARD	6850573
9	PCB SENSOR FOR O2 SENSOR	6850638
10	TEMPERATURE SENSOR, CPL.	8200486
12	MOISTURE SENSOR	8290393
13	PROTECTING CAP, without illustration	8290392
14	ANGLE,CPL.	8200497
17	ADHESIVE LABEL, without illustration	8303296
18	CABLE CLIP (2,4X92), without illustration	8712007
19	SHRINK HOSE PAC EX, without illustration	1202774
21	COVER,CPL., without illustration	8200597
22	CAPSULE FOR O2-DETECTOR (DW)	6850645

**BORDERING, CPL.****Picture 45**


Datei/File: 6141232.040


Item No.	Description	Order-Code
1	PROFILE	2M20195
3	WASHER A 4,3 DIN 125-A2	1327542
4	LOCKING PIN	2M20351
5	SCREW DWN 562-4X20-C-ST/A2P	1338382
7	FLAT PLUG DIN 46244-A6,3	1814079
8	SPLIT WASHER B4 DIN 127	1331345
9	HEXAGON NUT M4 DIN 934-A4/051	1328956
10	ADH. TAPE 12X1,6 PVC-FOAM SW X	1202421
11	LABEL, without illustration	2M20564
12	BUFFER	2M21887


**PLUG-IN UNIT 8000 IC.****Picture 46**


**PLUG-IN UNIT 8000 IC****Picture 47**



**PLUG-IN UNIT 8000 IC****Picture 48**


**PLUG-IN UNIT 8000 IC****Picture 49**

**PLUG-IN UNIT 8000 IC****Picture 50**

**PLUG-IN UNIT 8000 IC****Picture 51**


PLUG-IN UNIT 8000 IC

Picture 52




Item No.	Description	Order-Code
1-41	PLUG-IN UNIT 8000 IC with Skin Temp. Control, for SW 21.n only, without EPROM. Please order in addition sticker " $\pm 1.0\text{ }^{\circ}\text{C}$ " with the order no. 2M 21 921 for US versions in the USA and Canada.	2M21680
1-41	SERVICE EXCHANGE PLUG-IN UNIT (2M21680), for SW 21.n only, without EPROM Please order in addition sticker " $\pm 1.0\text{ }^{\circ}\text{C}$ " with the order no. 2M 21 921 for US versions in the USA and Canada.	2M21704
1-41	PLUG-IN UNIT 8000 IC with Skin Temp. Control, for SW 20.n only, without EPROM	8290443
1-41	SERVICE EXCHANGE PLUG-IN UNIT (8290443), for SW 20.n only, without EPROM	8290605
2	PCB ANALOG, <a href="#">SEE PICTURE 56</a>	
3	PCB DISPLAY/DISPLAY SKIN, <a href="#">SEE PICTURE 57</a>	
4a	KEY PAD, CPL., 8000 IC with Skin Temp. Control, for SW 21.0 and 20.n Please order in addition sticker " $\pm 1.0\text{ }^{\circ}\text{C}$ " with the order no. 2M 21 921 for US versions with SW 21.n in the USA and Canada.	2M21681
4b	STICKER " $\pm 1.0\text{ }^{\circ}\text{C}$ " FOR THEMOMONITORING	2M21921
4c	KEY PAD,CPL.(8000 IC with Skin Temp. Control), for SW 20.n only	8290467
5	ADH.TAPE 19X2-TESAMOLL 4700 WS	1209582
6	SPLIT WASHER B3 DIN127-X12CRNI	1309749
7	HEXAGON NUT M3 DIN934-M A4	1334913
8	NUT BM3 DIN 439-A4/051	1334743
9	BUTTON	8200619
15	FIXING BAR	8200899
17	BELLOW	8200927
19	FIBRE RING	D04766
23	FLAT PLUG DIN 46244-A6,3	1814079
25	SCREW	1287419
26	COUNTERS.SCREW AM3X6DIN963-A4	1287400
27	SCREW	1329987

Item No.	Description	Order-Code
28	CYLIND.HEAD SCREW AM4X20 DIN84	1328980
29	SCREW AM3X16 DIN 963-A2	1334557
30	WASHER A 4,3 DIN 125-A2	1327542
32	CH.HEAD SCR. AM3XDIN84-A4/051	1330799
33	PCB DISPLAY AIR, <a href="#">SEE PICTURE 53</a>	
34	PCB MOTHERBOARD	8290471
34a	RELAY BIPO 24/2 (MOTHERBOARD PCB)	1812653
35/34b	D PIA 6821 1MHZ DIL 40	1823361
37	PCB SWITCH	8200880
38	PCB POWER PACK, <a href="#">SEE PICTURE 54</a>	
39	PCB CPU INCUBATOR, <a href="#">SEE PICTURE 55</a>	
40	CABLE HARNESS, VOLTAGE REGULAT.	8306259
40a	NREGF 5V TO220	1830724
41	PCB DISPLAY/DISPLAY SKIN, <a href="#">SEE PICTURE 57</a>	
45	CAP	8200615
46	NREGF 5V TO220	1830724
47	SHR.D.PLAST.TUB.3,2-1,6X0,5	1200186
48	PCB MOTHERBOARD	8290471
48a	RELAY BIPO 24/2 (MOTHERBOARD PCB)	1812653
48b	D PIA 6821 1MHZ DIL 40	1823361
53	CABLE HARNESS PLUG-IN UNIT TO AGGREGATE	8200759
54	PCB CONTROLLER (for RS232), <a href="#">SEE PICTURE 58</a>	

**PCB DISPLAY AIR****Picture 53**

Item No.	Description	Order-Code
	PCB DISPLAY AIR (all versions)	8290411
D1	D PROG. KEYBOARD INTERFACE	1820877




**PCB POWER PACK****Picture 54**

Item No.	Description	Order-Code
1	PCB POWER PACK (all versions)	8290461
2	ACOUSTIC SIGNAL GENERATOR	8301708
4	NC STORAGE BATTERY 9V	8301856
5	ADDITIONAL FAN (not necessary)	see TSB # 17

## PCB CPU INCUBATOR


Picture 55



Item No.	Description	Order-Code
	PCB CPU INCUBATOR (all versions)	8290571
D6	MEMORY DECODER, CPL.	8304920
D8	D-RAM 16K+BATT+CLOCK(MK48T02B)	1828142
D12	SOFTWARE INCU.8000 IC 21.02 (for units with ThermoMonitoring only)	2M22331
D12	SOFTWARE 8000 IC 20.04 (for units w/o ThermoMonitoring only)	2M22327
D13	TIMER 6840	1813900
X	XFL SHORT CIRCUIT SOCKET 2POL.	1814508


## PCB ANALOG

Picture 56




Item No.	Description	Order-Code
PCB Analog for Incubator 8000 IC with ThermoMonitoring (units with SW 21.n):		
	PCB ANALOG, THERMOMONITORING, for all versions with SW 21.02 (replaces P/N 8290680, see TSB # 20)	2M22382
PCB Analog for Incubator 8000 IC without ThermoMonitoring (units with SW 20.n):		
	PCB ANALOG with Skin Temp. Control, for all version with SW 20.04 (replaces P/N 8200920, see TSB # 20)	2M22381
K1,K2	RELAY 5V/0,5A (K1/K2 of Analog PCB P/N 2M22400 and P/N 2M22401)	1813714

PCB DISPLAY / DISPLAY SKIN


**Picture 57**

Item No.	Description	Order-Code
1	PCB DISPLAY SKIN, for SW 21.n only	8290691
1	PCB DISPLAY, for - for Skin Temp. Module with SW 20.n only - for O2 and Humidity module (all versions)	8290431
2	D PROG.KEYBOARD INTERFACE	1820877

**INSTALLATION SET RS 232**
**Picture 58**


Item No.	Description	Order-Code
	INSTALLATION SET IC8000 RS 232	8290607
	BABYLINK CABLE BL.8000, without illustration	8306488


### Picture 60



Item No.	Description	Order-Code
1-9A, 10-13	CABLE CONNECTION THERMOMONITORING (Incubator 8000 IC with SW 21.n only)	2M21733
1-4, 6-7 9B-13	CABLE CONNECTION SKIN TEMPERATURE (YELLOW ONLY)	2M21738
6	XFL 7POL BUGE	1817604
9A	LABEL	2M21767
9B	LABEL	2M21759
13	CABLE CLIP (2,4X92)	8712007

**SKINTEMPERATURE SENSOR**


**Picture 62**



Datei/File: 6141232.053




Item No.	Description	Order-Code
For Incubator 8000 IC with ThermoMonitoring:		
1	SENSOR-SKIN TEMP-YELLOW (5 PACK), NEW VERSION, FOR CENTRAL SKIN TEMPERATURE	MX11000
3	SENSOR-SKIN TEMP-WHITE (5 PACK), NEW VERSION, FOR PERIPHERAL TEMPERATURE	MX11001
	HYDRO TEMP PROBE COVERS (100 pack) without illustration	DHNICU-00
For Incubator 8000 IC without ThermoMonitoring:		
7	ADHESIVE RING, 250 PIECES, FOR P/N 8200757	6850021
9	INTERFACE ADAPTER CABLE / SENSOR MALLINC.	2M20736
2	TEMP SENSOR, SINGLE USE, BLUE, FOR CENTRAL TEMPERATURE (WITH ADHESIVE PAD), TO BE USED WITH ADAPTER CABLE P/N 2M20736	D5030170
	HYDRO TEMP PROBE (BLUE ANGLED CONNECTOR), CAN BE USED WITHOUT ADAPTER 2M20736	DHNICU-32
	HYDRO TEMP PROBE COVERS (100 pack) without illustration	DHNICU-00

**MODIFICATION SET INCUBATOR IC, US LANGUAGE****Picture 63**

Datei/File: 6141232.054


Item No.	Description	Order-Code
1-6,9-13	MODIFICATION SET INCUBATOR 8000 IC USA	2M21937
12	O2 CONNECTION (USA)	8409401

**VENTILATION HOSE HOLDER****Picture 65**


Datei/File: 6141232.055

Item No.	Description	Order-Code
1-9	VENTILATION HOSE HOLDER, new version	8411075
1,5,10,11	VENTILATION HOSE HOLDER, old version	8411075
1	SCREW	2M15332
7	CLAMP	8411076
8	DISC	1335359
9	SCREW BZ 2,9X13 DIN 7971-A2	1334948
10,11	SET OF TUBE CLIPS 2	8406138

### Picture 66



Item No.	Description	Order-Code
1-2	DEPOSIT TRAY, COMPL.	2M20320
3	GASKET	2M20635
5	FLAT HEAD SCREW AM 2X4 DIN84	1287664
6-10	SET OF LABELS USA	2M20580

**MODIFICATION KIT CORNER PIECE**
**Picture 69**


Item No.	Description	Order-Code
1-10	MODIFICATION KIT CORNER PIECE	2M20722
9	STIFFENING	2M20214